






COTTON – Naturally in all sizes.

In each daily activity COTTON is the ideal companion. For all those who value natural fibres and a natural look, COTTON is the perfect choice. It is available in all sizes and many classic colours. Three different COTTON product lines means a perfect fit for all.










Circumference in cm	x-small	x-small PLUS	small	small PLUS	medium	medium PLUS	large	large PLUS
cT	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120
cH	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130
cG	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79
cG Grip top	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76
cF	41-52	44-55	44-55	49-60	49-60	53-64	53-64	59-70
cE	33-40	37-44	37-44	39-46	39-46	41-48	41-48	45-52
cD	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48
cC	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49
cB ¹	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37
cB	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29
cY	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38
cA	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29



Length in cm					
short	–	short	56-64	short	–
normal	35-40	normal	65-74	normal	65-74
long	41-45	long	75-84	long	75-84
 lD		 lG		 lG	






c = circumference
ℓ = length

		COTTON			COTTON XTRA	COTTON MED
CCL		1	2	3	1+2+3	1+2+3
Models						
	A-D Calf	with soft top	● ^{1) 2)}	● ²⁾	● ^{1) 2)}	⊙
		with grip top		● ¹⁾		
	A-G Thigh	with grip top	● ^{1) 2)}	●	■	
		with waist attachment (left or right)				●
	A-T Panty		■	●	■	
	A-T Panty Materna	for pregnant women and larger sizes	■	●	■	
	A-T Panty for men	with fly (left, right or horizontal)		■		
	open toe		●	●	●	●
	closed toe		■	●	■	⊙
	foot size one size			⊙		⊙
	maxi size					⊙
Grip tops						
	knobbed grip top (KGT)		●			
Colours						
	BLACK		●		⊙	
	ANTHRACITE				⊙	
	MARINE				⊙	
	CARAMEL		●			
	NATURE		●		⊙	●

CCL = Compression class
● Series production and made-to-measure articles
⊙ Series production
■ Made-to-measure articles

1) only in NATURE
2) only open toe

Material composition:
52% polyamide, 34% elasthan, 14% cotton

-  Wash daily using a laundry bag
-  Do not use fabric softener
-  Do not place in dryer
-  Do not iron
-  Do not dry clean

COTTON

