



TRADITIONAL – Many models, subtle tone.





This medical compression stocking for highest requirements is the first choice in serious venous problems both for women and also for men. TRADITIONAL is available in many sizes and model versions, as a series production or made-to-measure option, in compression classes 2, 3 and 4.












c in cm	x-small*	small	small PLUS	medium	medium PLUS	large	large PLUS	x-large*	x-large PLUS*
cT		up to 80	up to 80	up to 90	up to 90	up to 100	up to 100		
cH		up to 110	up to 110	up to 120	up to 120	up to 130	up to 130		
cG	40-50	42-52	50-60	50-60	58-68	58-68	66-76	61-71	69-79
cF	39-48	41-50	45-54	48-57	53-62	55-64	61-70	58-67	64-73
cE	29-34	32-37	36-41	35-40	40-45	38-43	44-49	41-46	48-53
cD	26-32	28-34	32-38	31-37	36-42	34-40	40-46	38-43	44-49
cC	28-31	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
cB ¹	21-24	24-29	24-29	29-32	29-32	32-36	32-36	35-39	35-39
cB	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32
cY	25-28	28-32	28-32	32-35	32-35	35-38	35-38	38-41	38-41
cA	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32

*only selected models



Length in cm							
short	33-37	short	56-60	short	61-70	short	61-70
long	38-44	long	61-70	long	71-80	long	71-80
 lD	 lF	 lG	 lG				

c = circumference
ℓ = length

TRADITIONAL				
CCL		2	3	4
Models				
 A-D Calf	with soft top	●	●	● ¹⁾
 A-F Half-thigh	with soft top	●	●	● ¹⁾
	with grip top	●	●	■
 A-G Thigh	with soft top	●	●	● ¹⁾
	with grip top	●	●	■
	with waist attachment (left or right)	●	● ¹⁾	■
 A-T Panty		● ¹⁾	■	■
 A-T Panty Bodyform	with compressing panty section	■	■	■
 A-T Panty Materna	for pregnant women and larger sizes	● ¹⁾	■	■
 A-T Panty for men	with fly (left, right or horizontal)	● ¹⁾	■	■
 open toe			●	
 closed toe			■	
Grip tops				
knobbed grip top (KGT)			●	
Colour				
BEIGE			●	

CCL = Compression class


● Series production and made-to-measure articles


■ Made-to-measure articles


¹⁾ only in PLUS


Material composition:


70% polyamide, 30% elastodiene

 Wash daily using a laundry bag

 Do not use fabric softener

 Do not place in dryer

 Do not iron

 Do not dry clean

